CAREvent®ALS

Handheld Automatic Transport Resuscitator





CAREvent® ALS

The CAREvent ALS Handheld Resuscitator provides trained individuals with a safe and effective means of providing artificial ventilation during respiratory and/or cardiac arrest.

The **CAREvent** ALS Handheld Resuscitator is a pneumatically powered, time/volume cycled resuscitator with the added feature of a Manually Actuated, Automatic Ventilation Override Button (Manual Button). No electrical supply or batteries are required to operate the device.

The single slider control for the simultaneous adjustment of ventilation frequency and tidal volume has six, preset, automatic and manual override settings for a range of patients from 20 kg. body weight to large adults.

During "mask ventilated" CPR the first breath is provided automatically followed by a manually actuated breath. This inhibits the automatic cycling for 20 seconds, allowing the rescuer to provide the 30 chest compressions recommended by the current guidelines for resuscitation without the ventilator attempting to provide the next breath.

If the Manual Button is not depressed again the automatic cycling will re-start at the selected rate and volume without any intervention by the operator being

An OFF position is provided so that the CAREvent ALS can be left attached to the gas supply with the supply turned on ready for immediate use.

The patient is protected from barotrauma due to over inflation by the audible, Pressure Relief System. This warns the operator of airway obstruction or airway restrictions. The excess delivered volume is discharged to atmosphere.

The **CAREvent** ALS Handheld Resuscitator also offers a "Demand Breathing" feature. This allows breathing patients to spontaneously breathe through the resuscitator on 100% oxygen, at their own desired rate and volume.

In addition, their inspiratory efforts will cause the automatic cycling to cease. Should the patient stop breathing again, the ventilator will restart cycling automatically at the selected rate and volume, without any intervention by the operator being required.

NOTE: To ensure proper ventilation of the patient in a weakened respiratory state, the automatic cycling shut off will only activate if the patient's inspiratory efforts are adequate.

Hand-held, lightweight, portable, simple to use and extremely durable, the CAREvent ALS Handheld Resuscitator is specifically designed for the demands of the emergency resuscitation and rescue environment.

Accessories available for the CAREvent ALS Handheld Resuscitator include the Deluxe Transport Ventilator circuit with PEEP port and the Gurney Clip (supplied with each unit), which allows the device to be attached to the hand rail of the gurney.

SPECIFICATIONS

Tidal Volume: 150 - 600 ml **Breaths Per Minute:** 20 - 10 I:E Ratio: 1:2 Automatic Flow Rate Range: 9 - 18 L/min Manual Flow Rate Range: As per Vt/BPM **Control Setting** Manual Overide Delay Time: 20 (+/-15%) Seconds Demand Breathing Flow Rate: 0 - 120 L/min Demand Breathing Triggering Pressure: -6cm H₂O@60 L/min

Auto Shut Off Delay Time: 5 - 8 Seconds Input Pressure: 45 - 70 PSI (3.1 - 4.8 Bar) Audible Pressure Relief Valve: 60 cm H₂0 -18°C to +50°C Operating Temperature: $(0^{\circ}F \text{ to } 122^{\circ}F)$

 -40° C to $+60^{\circ}$ C (-40° F to 140° F) Storage Temperature:

Relative Humidity for Storage

and Operating Use: 15 to 95% Input Connection: 9/16" DISS Patient Valve Dead Space: 8 ml.

Dimensions (inches): 5.5 x 2.5 x 2.9 (approx)

(millimeters): 140 x 63 x 73 (approx) 16 oz. / 450 gms (approx) Weight:

Ordering information: 01CV3000 CAREvent ALS Handheld Resuscitator

c/w 6 Foot Supply Hose and Universal Resuscitation Facemask

O-TWO MEDICAL TECHNOLOGIES INC. "INNOVATION IN RESUSCITATION"

7575 Kimbel Street, Mississauga, Ontario, Canada L5S 1C8 Toll Free: (800) 387-3405 Tel: (905) 677-9410 Fax: (905) 677-2035

E-mail: resuscitation@otwo.com

Website: www.otwo.com





DICAL TECHNOLOGIES INC